



“Q” PUMP INSTRUCTIONS

CERAMPUMP®

1506002

CONGRATULATIONS! The FMI “Q” PUMP LINE is one of the finest metering pump lines in the world. It features the unique valveless, variable, rotating and reciprocating pumping design of all FMI pump products and more. More consistent accuracy, more volume per stroke, more motor drive options and more diverse flow control options including a dial indicator, 4-20 milliamp stroke rate controller and dispenser.

SAFETY INSTRUCTIONS

! Before using any Fluid Metering, Inc. product read the following safety instructions as well as specific product specifications and operating instructions.

! **Warning!** Fire, electrical shock or explosion may occur if used near combustibles, explosive atmosphere, corrosive air, wet environment or submerged in fluid.

- Turn off the electrical power before checking pump for any problems.
- Connect motor, speed controllers, or any other electrical devices based on Fluid Metering Inc. specifications. Any unauthorized work performed on the product by the purchaser or by third parties can impair product functionality and thereby relieves Fluid Metering, Inc. of all warranty claims or liability for any misuse that will cause damage to product and /or injury to the individual.
- Power cables and leads should not be bent, pulled or inserted by excessive force. Otherwise there is a threat of electrical shock or fire.
- Replace any in-line fuses only with fuse rating as specified by Fluid Metering, Inc.
- When pump/drive is under operation, never point discharge tubing into face or touch any rotating components of pump.
- In a power down thermal overload cut-in condition, unplug or turn off power to pump. Always allow a cool down period before restarting; otherwise, injury or damage may occur.
- For 30 seconds after power is removed from pump/drive: do not touch any output terminals. Electrical shock may occur because of residual voltage.

! **Caution!** Fire, electrical shock, injury and damage may occur if not used in accordance with Fluid Metering, Inc. specifications and operation instructions.

- Do not put wet fingers into power outlet of unit.
- Do not operate with wet hands
- Do not operate drive assemblies that require a hard mount (to be bolted down) unless they are mounted per Fluid Metering, Inc. specifications, if not injury may occur and/or damage to unit.
- Do not touch any rotating pump or motor components: injury may occur.
- Do not run pump dry, unless designed for that service.
- Running dry is harmful to the pump, and will cause excessive heating due to internal friction.
- Check pump rotation and inlet/outlet pump port orientation before connecting power to pump. If not injury may occur.
- When pulling out cords from outlets do not pull cord, grasp plug to prevent plug damage or electrical shock.
- Fluid Metering, Inc. Drive Motors become **HOT** and can cause a burn. **DO NOT TOUCH!**

INSTALLATION & OPERATING TIPS

1. CLEAN FLUIDS. Abrasives in the pumped fluid may damage cylinder and piston surfaces and should, therefore, be avoided. Carbon cylinder liners and stainless steel pistons are particularly susceptible to abrasion by particulate matter in the pumped stream. Ceramic piston/cylinder components are more tolerant of suspended solids except solids that tend to flock and impede the movement of the piston in the cylinder.

2. COMPATIBLE FLUIDS. Pump only fluids compatible with materials of construction of the pump head you have selected.

3. WET OPERATION. The pumped fluid provides surface cooling and lubrication to the piston and cylinder of your FMI PUMP. Therefore, avoid dry operation (except pumps specifically designated “gas pump”).

4. FLOW VOLUME AND DIRECTION. Angular deflection of the cylinder with respect to the zero point on the calibration scale of your FMI “Q” PUMP controls flow magnitude and direction e.g., with the cylinder pointer at 10 on the left scale, fluid will be passed from the right port to the left port at 100% of the maximum rated volume; with the pointer at 10 on the right scale, fluid will pass from the left port to the right port at maximum rate. Set at 5 on the scale, flow rate will be 50% of maximum; at 4, it will be 40%; at 3, 30%, etc., etc. The flow control setting may be changed (including flow reversal) at any time while the pump is operating or

idle. Slightly loosen the two thumb screws and turn the STROKE LENGTH ADJUSTMENT KNOB. Retighten thumb screws once the desired setting is reached.

5. PISTON SEALS. The R408 seals that keep your PUMP piston dry are not “just ordinary plastic discs.” They are precisely cut and hot formed from sheets of a chemically inert fluorocarbon, specifically formulated for resistance to wear, abrasion, heat and chemical attack.

Each R408 seal possesses an exceptional mechanical memory which allows it to maintain a relatively constant wiping pressure on the piston, compensating for seal wear as it occurs. Properly maintained in clean condition, the original seals on a FMI PUMP may be expected to last the life of the pump. If they are removed for any reason, they should be carefully cleansed of all foreign particles prior to re-assembly. Seal seats must also be free of particles. (please see para 18)

6. DIAL INDICATOR. (optional) The Dial Indicator Kit is for fine adjustment and continuous monitoring of your “Q” pump flow rate settings.

To adjust Dial Indicator equipped pumps:

1. Loosen thumb screws, turn STROKE LENGTH ADJUSTMENT KNOB, moving cylinder assembly to neutral (zero-flow position).

2. Adjust indicator pointers until they read zero on both dials.

3. You are now ready for fine setting by turning STROKE LENGTH ADJUSTMENT KNOB until you achieve desired flow rate on dial. To prevent system backlash always turn STROKE LENGTH ADJUSTMENT KNOB two turns or one full revolution of large dial beyond desired setting, then adjust back.

7. 4-20 mA CONTROL for automatic response to remotely generated 4-20 milliamp signals is standard on V200 controllers. The input can be either grounded or ungrounded. The current source connects to terminal posts mounted on front cover assembly of the STROKE RATE CONTROLLER. Be sure to observe correct polarity. For complete hook-up and operating information see page 10.

8. PRESSURE. Do not operate pump against head pressures in excess of design specification. Drive arm on piston may bend or break under overload and other irreparable damage may be suffered. **Check your fluid circuit before applying power to the pump!**

9. ELECTRICAL PROTECTION. All FMI PUMPS are positive displacement instruments and should be protected by lowest possible “slo blow” fuse or circuit breaker electrical arrangements. “QV” units come equipped with .75 amp fuses.

10. NOISE AT HIGH PUMP RATES. A metallic hammering noise during operation of your pump (particularly high speed units such as QB, QD, QDX, and QV) when pumping liquids indicates

presence of gas bubbles in the pumping chamber which are reducing pumping capacity and may be damaging cylinder walls. Such bubbles may be traced to 1) a poor seal at the suction fitting, 2) fluid vaporization (cavitation) or, 3) degassing of the fluid.

a) To correct suction fitting leaks in stainless steel pump heads, remove fitting and wrap two layers of Teflon tape (standard Lab plumbing variety, 1 to 2 mil thick x 1/2" wide) tightly into the threads of the fitting. Replace fitting in cylinder port, drawing threads tightly on the Teflon tape. (see para 16).

b) To eliminate vaporization and degassing noise, reduce suction load. This may be accomplished by: 1) Using the 3/8" dia. TUBE ADAPTER R412-2 supplied with each pump on the suction line of the pump head to increase inside diameter of the suction line (use 1/2" dia. TUBE ADAPTER R412-6K on -3 PHM's.); 2) reduction of suction lift height; 3) pressurization of suction supply container; 4) locating pump below supply source to permit gravity flow aid; 5) reduce viscosity of fluid by heating or thinning; 6) reduce flow rate by adjusting pump to lower setting on flow scale; 7) install FMI PD-HF PULSE SUPPRESSORS in suction and discharge lines.

Improvements in noise abatement and pump life can be gained by putting pulse suppression hardware in the plumbing circuits adjacent to the pump suction and discharge ports - particularly with high speed pumps that are plumbed with rigid tubing. Theory holds that if part of a generated pulse is resiliently stored, the part not stored is smaller and thus easier to get into motion; the stored part of the pulse dissipating behind the part that is in motion sustains motion, causing an undulating flow to be transmitted rather than a series of pulses. Result: less noise, less energy used and less agitation of the pumped fluid. So for pulse noise and vibration problems, put a little resilience in your circuit. There are a number of rather easy ways to do it:

c) The simplest method is to use resilient tubing between the pump and the fluid circuit. Experiment a bit with standard elastomers - viton, hypalon, gum rubber, soft vinyl or other. Use only unreinforced tubing (reinforcement takes away the resilience). **Always shield this type of arrangement** so that a possible tube rupture will not endanger people or equipment.

d) Another popular pulse suppression arrangement involves a gas bubble trap as described in the final sentences of para. 12. A bubble in such a vertical trap will suppress pulse shock and noise temporarily. However, since gas and a liquid in contact under agitated conditions seldom stabilize, the trapped gas may absorb into the passing liquid and disappear leaving no pulse suppression or the fluid may contribute to the gas quantity, overload the trap and cause random pumping errors as occasional bubbles enter the flow stream. This can be overcome by fitting a soft slug of closed-cell-plastic foam or a soft pillow of thin-wall plastic tubing (ends sealed) into the vertical dead end extension of the fluid line. The gas trapped in the foam or pillow will provide the required resilience but will not be absorbed by the flow stream.

e) Since each fluid and circuit exhibits differing characteristics, a bit of experimentation may be necessary. The results are usually worth the effort.

11. FOR BEST LOW FLOW PUMPING RESULTS: Use a pump having a maximum flow rating as near

to the desired flow rate as possible and keep suction and discharge pressures essentially constant (see para 13). FMI pumps using R479 Low Flow Kits or designated LF are specifically designed for low flow/low dead volume, 1/4-28 flat bottom fittings.

12. LOW FLOW BUBBLE PROBLEMS. A common cause of trouble in metering pump applications requiring low flow rates - a few milliliters per minute or less - is the seemingly inevitable gas bubble trapped in the pumping head of the metering pump. It expands on the suction stroke and contracts on the discharge stroke, allowing little, if any, liquid to pass through the pump. Such bubbles, though often attributed to leaks in pump seals, can usually be traced to gases released by the pumped fluid in response to pumping agitation or pressure/temperature changes. When so identified, this potential source of metering pump error can be effectively controlled in most fluid circuits.

The familiar bubbles that form on the inside walls of a tumbler of tap water after it stands for a period of time at room temperature demonstrate the typical liquid degassing that results from pressure reduction (water line pressure to atmospheric) and/or temperature elevation (from ground ambient to air ambient). In this case, the bubbles contain air, hydrogen, carbon dioxide or other gaseous materials carried in the water; only small quantities of vaporized water are present. Some liquids respond to agitation and/or pressure/temperature changes by chemically separating into liquid and gas fractions; others simply vaporize, physically changing from liquid to gaseous form. Examples of liquids releasing gas or changing from liquid to gaseous form in response to agitation and temperature/pressure changes are numerous in the modern technical environment and many techniques have been devised to compensate for or correct their presence.

The most common practices for bubble control employ:

a) pressure on the suction side of the pump circuit to encourage gas retention in the liquid or,

b) employ natural buoyancy of the bubbles to carry them away from or through the pump head.

To apply pressure on the suction side of the pump, locate the pump physically below the supply vessel. Each two feet of elevation difference represents pressure of approximately one pound per square inch (psi). Bubbles that do occur will return to the supply vessel by buoyant lift. This is called a positive suction or flooded suction arrangement. If it is necessary to draw liquid up from the supply vessel to the pump head, negative suction pressure must be contemplated - again, approximately 1 psi per two feet of lift. Most liquids will release some gas when held at negative pressure and since the volume of gas released is generally proportionate to the volume of liquid subjected to the negative pressure, suction line diameter should be kept small for small flows (except heavy, viscous or tacky liquids which require large flow area for mobility). A vertical dead-end extension of the suction line can be provided above the pump suction port to trap line-generated bubbles before they enter the pump. This extension should be liquid filled at the start of a pumping period. Stand the pump vertically by loosening the screws and repositioning the Multi-Position Tilt Stand Q650 so that pump is in standing position, or hang the pump vertically by its base key slots. The discharge port should now be above the suction port allowing bubbles that enter the pump head to pass directly through with buoyant assist. Discharge lines should be inclined upward from pump head and bubble traps should be purged as often as necessary to

assure liquid flow continuity.

13. SYSTEM ACCURACY FACTORS. Several interrelated factors are involved in the exceptional operating accuracy possible in systems using FMI PUMPS. Of primary concern are the following:

a) FMI PUMP DISPLACEMENT precision is based on a simplified positive stroke mechanism which has no secondary linkages to produce stroke to stroke mechanical errors and has no gravity actuated or spring loaded valves to introduce random valve seating errors. The single mechanical linkage component between the PUMP piston and its drive elements is a precision spherical bearing which transforms circular drive motion into elliptical thrust motion (reciprocation). The total mechanical clearance of this linkage is less than 0.1% of the maximum pump stroke length or, approximately 0.0003". Thus it may be said that PUMP displacement precision (stroke to stroke) is in the order of the mechanical linkage clearance; that is to say, stroke to stroke displacement is reproducible to less than 0.5% within the rated capacity of a given pump model.

b) FMI PUMP VALVING is performed by a flat in the piston which is mechanically aligned with one cylinder port during the suction portion of each stroke and with the other cylinder port during the discharge portion of each stroke. The flat alignment is controlled by the single drive bearing discussed in the preceding sentences. The valve action is therefore mechanically precise, and free of random closure variations.

c) FLUID SLIP, a term commonly used to describe the migration of fluid around the internal moving parts of gear, lobe and vane pumps, is the volumetric difference between physical component displacement and fluid through-put of a pump system. In the FMI PUMP, slip loss refers to the fluid which passes through the clearance space (approx. .0002") between the piston and the cylinder wall. Since this clearance represents a restrictive passage of essentially constant dimension, it will be readily seen that the slip rate is determined by viscosity, pressure and time: e.g. assuming constant fluid viscosity and pressure, slip will be a smaller factor in a high repetition rate pump (short time per stroke) than in a low repetition rate pump. As viscosity increases and pressure decreases, time (or repetition rate) becomes less a significant contributor to slip loss.

d) STROKE REPETITION RATE is directly related to drive motor speed which in turn is influenced by work load and electrical supply voltage, i.e., motor speed decreases when work load increases and when electrical supply voltage (115 Volts AC) decreases. This motor speed variation may amount to as much as 15% for work load variations between zero discharge pressure and maximum rated discharge pressure. A 10% voltage drop may result in as much as 20% motor speed reduction when the pump is operating against a significant head pressure.


e) THE FLOW STABILITY (precision) of an FMI PUMP is therefore principally related to consistency in fluid slip rate and stroke repetition rate and these functions in turn are related to external system load factors such as viscosity, differential pressure and electric line voltage; i.e., when load factors remain essentially constant, slip rate and repetition rate remain essentially constant; when viscosity increases, fluid slip rate and stroke repetition rate both decrease; when differential pressure increases fluid slip rate increases and stroke repetition rate decreases.

In short, FMI PUMP PRECISION is influenced by



IMPORTANT

RECOMMENDED FMI PUMP MOUNTING FOR MAXIMUM PERFORMANCE



For maximum pump performance, mount the pump with motor at 12 o'clock and pump head at 6 o'clock position. This orientation will allow air bubbles that enter the pumping chamber to directly exit thru buoyant assist. Discharge lines should be inclined upward from pump head.

fluctuations of fluid differential pressures, fluid viscosity and electric line voltage. When these factors are controlled predictably reproducible pumping precision better than 0.5% may be expected.



MAINTENANCE & REPAIR INSTRUCTIONS

14. LUBRICATION. Pump drive elements are provided with oil fittings at all appropriate points. Use high grade machine oil at regular intervals. A dab of good grease on PISTON DRIVE PIN just before it is inserted into the RADIAL BEARING in the SPINDLE ASSEMBLY does a world of good for the bearing and pin.

15. CHANGING THE LIQUID END OF YOUR FMI "Q" PUMP. (figures 1,2).

- a) To remove "Q" Pump Head Module (QPHM):
1. Turn power off.
 2. Rotate STROKE LENGTH ADJUSTMENT KNOB to position PUMP HEAD CARRIER Q410-2 all the way to the extreme right or left of scale (fig. 1).
 3. Rotate SPINDLE ASSEMBLY Q424 to place PISTON DRIVE PIN at 3 or 9 o'clock position (facing cylinder head).
 4. Loosen two KNURLED NUTS.
 5. Lift QPHM and draw gently away from Q424 SPINDLE ASSEMBLY (fig. 2).
 6. Move QPHM up and to the left while slipping DRIVE PIN out of RADIAL BEARING.

- b) To replace QPHM:
1. With PISTON ASSEMBLY R423 extending approximately 1-1/4" from CYLINDER NUT R406K and DRIVE PIN in the 3 or 9 o'clock position, insert PIN into the RADIAL BEARING in the SPINDLE ASSEMBLY.
 2. Slide QPHM into position on BASE making sure locator on bottom of CARRIER Q410-2 drops into slotted portion of FOLLOWER Q454-2 on base.
 3. Tighten KNURLED NUTS on Q616 assembly.





16. CHANGING FITTINGS ON STAINLESS STEEL PUMP HEADS IMPORTANT!

Fittings screwed too tightly into stainless steel pump heads will contact port seals and may cause piston/cylinder damage. Use extra layers of Teflon tape on threads when necessary to avoid such excessive penetration.

17. CLEANING PUMP HEAD. Routine flushing with solvent before shut-down will suffice for most applications - set pump for maximum stroke and operate until solvent appears clear at discharge port. If periodic teardown for detail cleaning is required, remove parts with care to avoid damage to piston, cylinder and gland. Wipe all parts with lintless oil saturated cloth. Operate by hand after reassembly to assure free movement of parts prior to application of power.

17.1 CAUTION! Ceramic piston/cylinder sets are particularly sensitive to neglect and may "freeze" if allowed to dry out without adequate cleansing. Some users actually remove the piston from the cylinder after solvent cleaning and store the component parts in disassembled condition until the pump is again required. Others fill a loop of flexible tubing with fluid that will thin or neutralize the last fluid pumped. They then connect one end of the tube to the pump suction port, the other to the discharge port. With this loop positioned above the pump head, the ceramic surfaces and seal areas will stay moist and mobile for extended idle periods. If, however, a piston (ceramic or stainless steel) does freeze in the cylinder, **DO NOT TRY TO FORCE IT FREE!** Be gentle. Try to remove the pump head (refer para 15) from the base assembly so that the whole assembly can be soaked in a suitable solvent. If the head is not conveniently




removable, the tube loop discussed in the prior paragraph may permit solvent to dissolve the "frozen" residue in reasonable time. If all else fails, pack it all up and ship it back to the factory with a note telling us what you think might be left inside - sometimes we're lucky. Having a spare pump head on hand in case of emergency is always a good idea.

17.2 SANITARY SERVICE. FMI sanitary pump heads, designated SAN are designed to conform with the cleansing standards of the U.S. Food and Drug Administration.

- To clean individual component parts:
- a) Unscrew 2 CYLINDER NUTS R406-S and 2 PORT NUTS 110348-P.
 - b) Remove piston and seal components released by step a) above;
 - c) Use fingers to wiggle LINER R407-C slipping it from CYLINDER CASE R405-SAN;
 - d) Cleanse and sterilize component parts as per government regulations for parts manufactured of Alumina Ceramic, Type 316 Stainless Steel and Teflon. After cleaning, reassemble the components as follows:
 - e) Assemble piston/seal components as per instructions para 18 a,b,c;
 - f) Orient CYLINDER LINER R407-C in CYLINDER CASE R405-SAN in such a manner that the side port flats on CYLINDER LINER R407-C should be visible through the side port openings on CYLINDER CASE R405-SAN.
 - g) Place one each port seal component R412-T into each side port of CYLINDER CASE R405-SAN using finger pressure across the two ports to assure flush seating of the seal surfaces on the cylinder liner flats.
 - h) Apply PORT NUT 110348-P to each port, finger tight.
 - i) Carefully insert piston into cylinder liner until threads of CYLINDER NUT R406-S meet threads of CYLINDER CASE R405-SAN. Tighten NUT R406-S finger tight.

Continued on Page 12

FLUOROCARBON PUMP HEAD MODULE



* REPLACEMENT CYLINDER CASE ONLY SOLD WITH CYLINDER LINER

IMPORTANT!


WHEN ORDERING REPLACEMENT PARTS PLEASE MENTION MODEL AND SERIAL NUMBERS OF THE PUMP ON WHICH THEY WILL BE USED.

REPLACEMENT CYLINDER CASE (R405-K) ONLY SOLD WITH CYLINDER LINER

CERAMIC PISTONS & CYLINDER LINERS SOLD AS SETS ONLY


PLEASE NOTE: FMI PROVIDES A 48 HOUR TURN AROUND REPAIR SERVICE (SEE PAGE 7 FOR DETAILS)

STAINLESS STEEL PUMP HEAD MODULE



CERAMIC PISTONS & CYLINDER LINERS SOLD AS SETS ONLY

SANITARY PUMP HEAD MODULE



PUMP HEAD INFORMATION

THINGS YOU SHOULD KNOW BEFORE ORDERING THE FLUID HANDLING PORTION (or parts thereof) OF AN FMI MODEL Q PUMP, i.e. the removable assembly that is often referred to as the **liquid end**, the **piston/cylinder set**, or the **Pump Head Module (PHM) assembly**.

To accommodate the diverse demands of laboratory and industrial pumping applications, FMI's Q-Pumps are ar-

ranged to accept interchangeable PHM's with differing chemical and mechanical characteristics. There are a number of these assemblies from which the pump user may select in solving difficult fluid pumping problems. To simplify the selection process, each individual type of fluid handling assembly (piston/cylinder set) is designated by an alpha-numeric code which permits direct identification of pertinent factors.

Q-PUMP DESIGNATIONS

Example:

Q-1-SAN denotes a sanitary construction pump head assembly having a 1/4" diameter ceramic piston, ceramic cylinder liner and a 316 stainless steel cylinder case.


PLEASE NOTE:

- When the last three of the above shown categories are **not** called out, standard materials and construction are implied (CKC).
- If **SAN** is designated in the last category, the preceding three categories will be designated CSC as required by FDA specifications.

3. Piston seals made of an inert material, Rulon AR, are used in all FMI type Q Pump Heads except high speed pumps fitted with stainless steel pistons (these use Rulon J seals) and sanitary (SAN) pump heads which require virgin Teflon (T) seals. Specify AR, J or T after the pump head code designation to order pump heads assembled with other than standard seals.

4. Ceramic cylinder liners should be used with ceramic piston only.

Pump Head Material Configuration



PUMP HEAD MODULE CODES									
PISTON SIZE CODE	MATERIALS OF CONSTRUCTION								
	CKC	CKY	CSC	CSY	SAN	SKY	SSY	STY	CTC
RH00		N/A	N/A	N/A	N/A				
RH0									
RH1		N/A	N/A	N/A	N/A	N/A	N/A	N/A	
RH Pump Heads: 1) Require an RH/Q kit for use with Q drives									
Q0	N/A	N/A	N/A	N/A	N/A			N/A	N/A
Q1								N/A	N/A
Q2									
Q3		N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
WETTED PARTS	CERAMIC PVDF	CERAMIC PVDF CARBON	CERAMIC 316 SS	CERAMIC 316 SS CARBON	CERAMIC TEFLON	316 SS PVDF CARBON	316 SS CARBON	316 SS TEFLON CARBON	CERAMIC TEFLON
MAX. TEMP	212 F	212 F	350 F	350 F	350 F	140 F	140 F	140 F	212 F
OPTIONS (add Option Code & cost to Pump Module for complete price and part number)									
LF	N/C	N/C	N/A	N/A	N/A	N/C	N/A	N/C	N/C
W		N/A		N/A	N/A	N/A	N/A	N/A	N/A
WT	N/A	N/A		N/A	N/A	N/A	N/A	N/A	N/A
TC	N/A	N/A		N/A	N/A	N/A	N/A	N/A	N/A
R479	N/A	N/A			N/A	N/A		N/A	N/A
S("Q" Only)								N/A	N/A

See Materials of Construction section in the Fluid Metering catalog for more information on wetted parts *316 Stainless Steel cylinder cases accept 1/4 NPT male fittings.

When ordering specify: Piston Size Code + Material of Construction + Option Code (Q1 + CKC + W = Q1CKC)

Q PUMP HEAD MODULE - PARTS PRICE LIST

PART NO.	DESCRIPTION	PART NO.	DESCRIPTION
Q410-2	CARRIER	R412-1	ADAPTER, S.S. 1/4 NPT TO 1/4" I.D. TUBE
Q410-3	CARRIER STAINLESS STEEL	R412-2	ADAPTER, S.S. 1/4 NPT TO 3/8" I.D. TUBE
Q665	SPACER	R412-0K	ADAPTER, FLUOROCARBON 1/8" I.D. TUBE
R405-K	CYL CASE, FLUOROCARBON, DBL CAP	R412-1K	ADAPTER, FLUOROCARBON 1/4" I.D. TUBE
R405-1K	CYL CASE, FLUOROCARBON, STANDARD	R412-2K	ADAPTER, FLUOROCARBON 3/8" I.D. TUBE
R405-S	CYL CASE, STAINLESS STEEL	R412-5K	ADAPTER, FLUOROCARBON 1/4-28 THD
R405-1S	CYL CASE, SS, TEMP GLAND	R412-6K	ADAPTER, FLUOROCARBON 1/2" I.D. TUBE
R405-SAN	CYL CASE, SANITARY STAINLESS STEEL	R412-07	ADAPTER, S.S. 1/4" I.D. TUBE, SAN
R405-SG	CYL CASE, SS ISOLATION GLAND	R412-08	ADAPTER, S.S. 3/8" I.D. TUBE, SAN
R406-K	GLAND NUT, FLUOROCARBON	R412-1T	ADAPTER, TEFLON 1/4" I.D. TUBE, SAN
R406-0K	GLAND NUT, FLUOROCARBON 1/8" PISTON	R412-2T	ADAPTER, TEFLON 3/8" I.D. TUBE, SAN
R406-1K	GLAND NUT, FLUOROCARBON 1/4" PISTON	R413	CYLINDER HEAD SEAL
R406-2K	GLAND NUT, FLUOROCARBON 3/8" PISTON	R413-1	CYLINDER HEAD SEAL, R405-1K
R406-3K	GLAND NUT FLUOROCARBON	R419	CYLINDER PORT SEAL, R405-S
R406-S	CYLINDER NUT R405-S, R405-SAN	R423-1C*	PISTON
R407-0	CYLINDER LINER	R423-2C*	PISTON
R407-1	CYLINDER LINER	R423-3C*	PISTON
R407-2	CYLINDER LINER	R423-0S	PISTON
R407-1C*	CYLINDER LINER	R423-1S	PISTON
R407-2C*	CYLINDER LINER	R423-2S	PISTON
R407-3C*	CYLINDER LINER	R479	LOW FLOW KIT (FITS R405-S)
R408-1A	LIP SEAL, RULON AR 1/4"	R489	SPACER, SANITARY CYLINDER GROUP
R408-2A	LIP SEAL, RULON AR 3/8"	R490-1	PORT SLEEVE, SANITARY CYLINDER GROUP
R408-0J	LIP SEAL, RULON J 1/8"	R520-1T	ADAPTER, TEFLON 1/4" O.D. STRAIGHT, SAN
R408-1J	LIP SEAL, RULON J 1/4"	R520-2T	ADAPTER, TEFLON 3/8" O.D. STRAIGHT, SAN
R408-2J	LIP SEAL, RULON J 3/8"	H476-K	SML TUBE ADAPTER SET, 1/8" O.D. (OPTL)
R408-1T	LIP SEAL, TEFLON 1/4"	110017-10	SCREW, SKT HD CAP 1/4-20 X 5/8"
R408-2T	LIP SEAL, TEFLON 3/8"	110311	NUT, SQUARE 1/4-20
R408-3A	LIP SEAL, RULON AR 1/2"	110348-P	SAN PORT NUT, PLASTIC
R408-3T	LIP SEAL, TEFLON 1/2"	110348-S	SAN PORT NUT, STAINLESS STEEL
R409-0	GLAND WASHER, TEFLON 1/8"	110384-K	FERRULE NUT 1/4, KYNAR
R409-1	GLAND WASHER, TEFLON 1/4"	110384-T	FERRULE NUT 1/4", TFE
R409-2	GLAND WASHER, TEFLON 3/8"	110847-01	BARBED FITTING 10-32 UNC-2B 1/8" ID TUBING
R409-3	GLAND WASHER TEFLON 1/2"	200182-01	CYL CASE, SS HIGH TEMP GLAND



* CERAMIC PISTONS & CYLINDER LINERS SOLD AS SETS ONLY

PUMP HEAD REPAIR SERVICE

Should your FMI PUMP or PUMP SUBASSEMBLY need repairs, you should first contact the Customer Service Department for a **SERVICE SAFETY ASSURANCE FORM** and a **RETURN AUTHORIZATION NUMBER**, you may then ship it to us post-paid. You will be charged a flat \$29.00 service fee for each Pump Head Module returned plus the cost of parts used and the postage required to return it to you. The Pump Head Module service fee includes inspection, clean and test to agreed specifications and the replacement of Lip Seals and Gland Washer.

If after examining a returned pump or subassembly, the FMI Service Department estimates the total cost to be in excess of \$60.00, you will be requested to approve the charges before repair work is started.

PARTS ORDERS

MINIMUM ORDER \$15.00 DOMESTIC or FOREIGN (Invoice price exclusive of shipping)

SHIPPING

Parts and repair orders will be shipped via United Parcel Service or U.S. Postal Service unless other means are specified.


ALL PRICES ARE QUOTED IN U.S. DOLLARS, FOB SYOSSET, NY - Subject to change without notice.

For Additional Information Call - Toll Free **800-223-3388** or email us at: **pumps@fmipump.com**




FMI MODEL Q PUMPS


PARTS IDENTIFICATION SHEET IN-Q431-04




QD PUMP DRIVE MODULE



QG PUMP DRIVE MODULE




QSY PUMP DRIVE MODULE



QB, QV, QVG PUMP DRIVE MODULE

Q421-3 BASE ASS'Y, QB, QD, QV
Q421-4 BASE ASS'Y, QSY
Q421-6 BASE ASS'Y, QG, QVG50



ACCESSORY






FMI MODEL Q PUMPS


PARTS IDENTIFICATION SHEET IN-Q431-04



Q421-5 BASE ASS'Y, QSYX
Q421-7 BASE ASS'Y, QDX



QP PUMP PEDESTAL



IMPORTANT!

WHEN ORDERING REPLACEMENT PARTS, PLEASE MENTION MODEL & SERIAL NUMBERS OF PUMP ON WHICH THEY WILL BE USED.

PLEASE NOTE FMI PROVIDES A 48 HOUR TURN AROUND REPAIR SERVICE (see page 11 for details)

FOR ADDITIONAL INFORMATION
 CALL TOLL FREE 800-223-3388 or
 516-922-6050 ● FAX 516-624-8261
 email: pumps@fmipump.com



V200 STROKE RATE CONTROLLER

INSTALLATION & OPERATING TIPS

1. CONTROL FUNCTION: The FMI stroke rate controller, V200 used with Models QV, Q2V, QVG50, and RHV pumps has four switches to control its functions: a) A two position power on-off switch. b) A two position direction switch. c) A Digital Potentiometer for control of percent of flow. A setting of 99.9 yields approximately 100% of maximum flow. d) A two position control method switch. In the up or manual position it allows for direct control of percent of flow by setting the Digital Potentiometer. In the down or 4-20 mA position, control is switched to the 4-20 mA terminal strip inputs and the Digital Potentiometer is disabled. Stroke rate is then controlled via a customer supplied DC current source with an input impedance of approximately 500 Ohms.

2. FUSING: The motor and its control are protected against overloads by a current limit circuit and additional protection is provided by a .75 amp fuse on input line (hot side). Replacement fuses must be .75 amp Slo Blow.

3. TRIM POT ADJUSTMENTS: The V200 trim pots are factory set and should not be adjusted. 0 = minimum and 1800 = maximum strokes per minute. If other than standard settings are desired please note Trim Pot Adjustment Methods.

4. TRIM POT ADJUSTMENT METHODS: Before adjusting control unit, carefully check that all con-

nections are correct. Upon removal of V200 cover, use a non-metallic screwdriver to avoid shorting of the control circuit. After adjustment and before replacing the cover, be sure no foreign materials are on the circuit board.

A) Minimum Speed Adjustment- The MIN. SPD adjustment sets the speed that the motor will attain when the Digital Potentiometer is set to 0% of flow which is the factory preset. To adjust minimum speed, set the Digital Potentiometer to 0% flow (00.0) and advance the MIN. SPD adjustment until the desired speed is achieved. The range is 0% to approximately 50% of the control rated output voltage (90 VDC). This adjustment should be made before setting the MAX. SPD adjustment. There is some interaction between settings.


B) Maximum Speed Adjustment - This sets the speed that the motor will attain when the Digital Potentiometer is set to 99.9% flow (99.9). Maximum speed is set at the factory to 1800 spm. The setting of the MAX. SPD adjustment has no effect when operating in the SIGNAL mode. To set a different maximum speed, turn the Digital Potentiometer to 99.9% of flow and adjust the MAX. SPD adjustment. The maximum speed may be set as low as 50% of the control rated output voltage (90 VDC).

C) Input Signal Scaling- Two adjustments may be made to calibrate the V200 to follow its signal source. First, apply the 20 mA signal into the control and adjust the SIG ADJ. trimpot to bring the motor to 1800 spm. Then apply a 4 mA signal to the control and adjust the LINEARITY trim pot to bring the motor to the speed desired at that low signal. In SIGNAL mode, verify that the input signal is at 4 mA. In MANUAL mode, verify that the Digital Potentiometer is set to 0% of flow (00.0). Slowly increase the signal input, or slowly increase the percent of flow setting. The motor must start slowly and increase its speed in approximate proportion to the increasing speed setting. If the motor rotates opposite to the intended direction, shut down the control and reverse the connections made to terminals A1 and A2, located on the edge of the circuit board nearest the digital potentiometer. Restore power to the control.


D) IR Comp Adjustment- This controls the degree to which the V200 compensates for changes in motor load to maintain essentially constant motor speed over the load range. It is factory set with the rated horsepower motor.

5. POWER REQUIREMENTS: The V200 stroke rate controller requires a 115 VAC, 50/60 Hz power source in order to function. For usage at 230 VAC @ 50 Hz FMI recommends the use of a step-down transformer rated at 150 VA minimum.

MANUAL CONTROL RESPONSE CURVES



INPUT SIGNAL RESPONSE CURVES



Q PUMP DRIVE MODULE - PARTS PRICE LIST

PART NO.	DESCRIPTION	PART NO.	DESCRIPTION
Q421-1	BASE ASS'Y QP	110290-10	SCREW, SKT HD #8-32 X 5/8"
Q421-3	BASE ASS'Y QB, QD,QV, QBG	110290-22	SCREW, #8-32 X 1 3/8"
Q421-4	BASE ASS'Y QSY	110297-6	DRIVE SCREW , 8-32 X 3/8"
Q421-5	BASE ASS'Y QSYX	110306-1	MOTOR, QSY 115V, 60Hz
Q421-6	BASE ASS'Y QG, QVG50	110306-2	MOTOR, QSY 230V, 50Hz
Q421-7	BASE ASS'Y QDX	110314	WASHER, #10 INT. LOCK
Q424	SPINDLE ASSEMBLY QP, QDX	110336-2	MOTOR 230 VAC 50Hz QSYX
Q447	BUSHING, FLOW CONTROL	110337-1	RESISTOR (110336-1 115V, 60Hz)
Q454-3	FOLLOWER	110337-2	RESISTOR (110336-2 230V, 50Hz)
Q464	STROKE LENGTH ADJUSTMENT KNOB	110338-1	CAPACITOR (110336-1 115V, 60Hz)
Q485	DIAL INDICATOR KIT	110338-2	CAPACITOR (110336-2 230V, 50Hz)
Q616	RETAINER ASSEMBLY	110355-1	MOTOR, 12 VDC QB
Q647	BASE MOUNTING KIT	110373	DRIVE COUPLING
Q662	TILT STAND KIT	110376-16	BOLT, HEX HD 3/8 -16 X 1"
		110417-4	STANDOFF, HEX 8-32 X 1/40 FEM.
R424-1	SPINDLE ASSEMBLY QD, QG, QB, QV, QVG	110454	CAP NUT 8-32
R424-6	SPINDLE ASSEMBLY QSY, QSYX	110471	LOCK NUT 1/4-20
R475	FOOT QSYX	110491	WASHER, FLAT 1/4" STD
		110505	WASHER, SHOULDER, NYLON
P56C	FACE MOTOR ADAPTER KIT	110506-8	SCREW, SKT HD CAP 10-32 X 1/2"
V109	MOTORASSEMBLY QV	110509	ADAPTER, PC56C
V109G	MOTORASSEMBLY QVG50	110522	MOTOR, X-PROOF QDX 115/230 VAC 60 HZ
		110542-6	SCREW, HEX HD NYLON 8-32 X 3/8"
110049	WASHER #8 INT. LOCK	110544	SPEED KNOB
110258	FOOT, RUBBER	110579	WASHER, FLAT 3/8"
110262	MOTOR, QD PUMP 115V, 60Hz	110580	WASHER, LOCK INT. 3/8"
110262-2	MOTOR, QD PUMP 230V , 50/60 Hz	110608	SPRING AND BRUSH SET, QVG50
110286-6A	MOTOR, QG 6 RPM 115V, 60Hz	110609	BRUSH CAP , QVG50
110286-20A	MOTOR, QG 20 RPM 115V, 60 Hz	110671-4	SCREW, SET 8-32 X 1/4"
110286-50A	MOTOR, QG 50 RPM 115V, 60 Hz	110837-01	SPRING BRUSH SET 12 VDC - GRAY MOTOR
110286-150A	MOTOR, QG 150 RPM 115V, 60 Hz	110838-01	SPRING BRUSH SET 90 VDC - GRAY MOTOR
110286-400A	MOTOR, QG 400 RPM 115V, 60 Hz	110839-01	BRUSH CAP 12, 90 VDC - GRAY MOTOR
110290-8	SCREW, SKT HD 8-32 X 1/2"		

Q PUMP DRIVE MODULE

PUMP DRIVE MODULE CODES												
PUMP DRIVE	QD	QG	QSY	QB	QBG	QV	QVG50	Q2V	Q2B	QP	QSYX2	QDX

PUMP DRIVE REPAIR SERVICE

Should your FMI PUMP DRIVE or SUBASSEMBLY need repairs, you should first contact the Customer Service Department for a **RETURN AUTHORIZATION NUMBER**; you may then ship it to us post-paid. You will be charged a flat \$21.00 service fee for each Pump Drive Module returned plus the cost of parts used and the postage required to return it to you.

If after examining a returned pump or subassembly, the FMI Service Department estimates the total cost to be in excess of \$60.00, you will be requested to approve the charges before repair work is started.

PARTS ORDERS

MINIMUM ORDER \$15.00 DOMESTIC or FOREIGN (Invoice price exclusive of shipping)

SHIPPING

Parts and repair orders will be shipped via United Parcel Service or U.S. Postal Service unless other means are specified.

ALL PRICES ARE QUOTED IN U.S. DOLLARS, FOB SYOSSET, NY - Subject to change without notice.

For Additional Information Call - Toll Free **800-223-3388** or email us at: pumps@fmipump.com

- j) Place head seal components R413 and R489 into the head end of CYLINDER CASE R405-SAN, apply CYLINDER NUT R406-S.
- k) Wrench tighten each PORT NUT 110348-P 1/6 turn.
- l) Wrench tighten each CYLINDER NUT R406-S 1/2 turn, piston end first.

18. PISTON SEAL REPLACEMENT (please see para 5.) When R408 SEALS are replaced, the following procedure should be followed: (please see figs. 3,4)

- a) Place GLAND NUT R406 and GLAND WASHER R409 on PISTON ASSEMBLY R423.
- b) First "form" lip of LIP SEAL around piston by gently placing a LIP SEAL R408 on piston, **lip side last**. Carefully rotate the seal on the piston

to avoid damage to the lip while passing over the flat to the piston neck. **Then remove seal and reverse lip direction** (fig. 3).

c) Gently place one "formed" LIP SEAL R408 on piston, **lip side first**, carefully rotating the seal on the piston until it is past the flat and on the piston neck.

d) Gently place one LIP SEAL R408 on piston, **lip side last**. Carefully rotate the seal on the piston to avoid damage to the lip while passing over the flat to the piston neck (fig. 4).

e) Insert piston into cylinder approximately one inch.

f) Apply GLAND NUT R406 to cylinder threads and tighten.

19. PISTON SEAL SETTING. After installing new

lip seals (part R408) in pump head it is recommended that the seals be set (formed in place) by fluid pressures generated by pump action. To accomplish this:

a) Operate the pump spindle clockwise for 10 or 20 strokes at maximum setting, handling water (left to right mode facing pump head) with suction line blocked or pinched off. This will create a vacuum in the pump head, permitting atmospheric pressure to shape the outer seal member tightly around the piston.

b) Reverse the pumping direction (pump head angle reversal) and intermittently block the line leading from the left hand port. This will generate pressure in the seal area of the pump head, causing the inner seals to form intimately around the piston.

IMPORTANT

WHEN ORDERING REPLACEMENT PARTS PLEASE MENTION MODEL AND SERIAL NUMBERS OF THE PUMP ON WHICH THEY WILL BE USED.

Record your Model and Serial Numbers below. This will be useful when ordering replacement parts.

PUMP HEAD MODULES (PHM)

PUMP DRIVE MODULE (PDM)

PHM MODEL NUMBER

PHM SERIAL NUMBER
QB _____

PDM MODEL NUMBER

PDM SERIAL NUMBER
QA _____

PLEASE NOTE:

Pistons and cylinders should be dimensionally mated at the factory. For most satisfactory performance, they should be ordered as mated sets.

NOTES:



FLUID METERING, INC.

5 AERIAL WAY, SUITE 500, SYOSSET, NY 11791
(800) 223-3388 • (516) 922-6050 • FAX (516) 624-8261 • <http://www.fmipump.com>

